

THE 100% NO TOUCH  
TITANIUM HOLDER

# THE HOLDER FOR A KLEAN DENTAL SURGERY



## PROTECTIVE:

The implant is held thanks to a bespoke titanium ring and its own covering/healing screw.



## NO TOUCH:

Surgeons / Scrub nurses can remove directly the implant from the holder without touching it, thus avoiding any contamination.



## SMART:

Compatible with the «Clean Implant» label. Suitable for right/left handed use.




## CUSTOMIZABLE:

SoKlean is available in 9 different shades, and its titanium ring can be customized with a wide range of colors. It can be packaged in pouches, blisters and SoBliss®.



## ADJUSTABLE:

SoKlean® can hold several implant sizes thanks to its two different versions and to the unidirectional slider.



SoKlean® and ring are available in different colors.



## Implant size

### SoKlean M

Dmax: 5,5 Lmax: 18  
Dmin: 3 Lmin: 6

### SoKlean L

Dmax: 7 Lmax: 24  
Dmin: 3 Lmin: 11


## Soklean® sterile packaging



## User guide



**TILT**  
LEFT OR RIGHT



**HOLD**  
THE BODY



**PICK**  
SAFELY



**PUSH**  
THE SLIDER



**GET**  
THE SCREW



www.selenium-medical.com | contact@selenium-medical.com

Selenium Medical: 40 rue Chef de Baie, 17000 L Rochelle FRANCE | Tel: +33 (0)5 46 44 40 82

Selenium America: 6201 Fairview Road, Suite 200, 28210, Charlotte NC USA | Tel: +1 980 215 8506



001.00.04 BR021.03.05.2019