



STERILE
SINGLE OR DOUBLE
BLISTER PACKAGING



THE REVOLUTIONARY BLISTER PACKAGING



EASILY INTEGRABLE IN YOUR CLEAN ROOM



VIDEO DEMO



Size L



WORLDSTAR
WINNER 2019



SoBliss® introduces the blister packaging V2.0. No changes have been implemented around the blister concept for decades. SoBliss opens a new era in the blister world. Thanks to the vertical all-around sealing and the hinge, a puncture proof Tyvek® system has been created. Medical Devices can now be packaged in a new rigid 360° transparent double blister.

SAFE

- » Single/Double/Triple sterile barrier
- » Validated gamma irradiation
- » Compatible ETO sterilization

OPTIMIZED

- » Various implant configurations possible
- » Vertical sealed double blister with hinge

BEST USE

- » Traceability improvement
- » Fast opening sleeves
- » Easy product identification

LIGHT

- » Smaller than a traditional blister
- » Shipping and storage cost saving

INNOVATIVE

- » Puncture-proof Tyvek® system
- » Worldwide patented*

SMART

- » Implant stability thanks to the SoTab® /SoEasy® holders
- » Single use warranty with tamper proof sleeve
- » Single use tyvek band opening

SURPRISING

- » 360° implant visibility
- » Ultra flat packaging

USABILITY STUDY (AVAILABLE UPON REQUEST)

Nurse Panel - AORN 2018 exhibition
SoBliss evaluated 8.8/10



so easy®

100% NO TOUCH
HOLDER FOR
SCREW

so bliss® Size M
FLAT



CLAMSHELL
CUSTOMIZED
HOLDER

so bliss® Size S
FLAT



ISO 11607
ISO EN 868-5




FR patent: 1354717/ 1660342
US patent: 14 892 059

Size S


Size M

Size L

EXTERNAL CAP
Protection & Identification



INTERNAL BLISTER
Tyvek band Sealing



FIT IN
⇓

EXTERNAL BLISTER
Tyvek band Sealing



ALL DIMENSIONS IN MILLIMETERS